

TAKE THE FLIGHT OF SUCCESS WITH VIGA  
GEMILANG

B-G-8, RAWANG IDAMAN 1,  
JALAN IDAMAN 1,  
TAMAN RAWANG IDAMAN, 48000 RAWANG SELANGOR

## CiTest Malaria Rapid Test

### Your True Companion in Malaria Screening

#### CITEST Malaria Rapid Tests Offer:

Malaria P.f. (PfHRP-2) Detection  
Malaria P.v. (Pv-pLDH) Detection  
Malaria Pan (Aldolase) Detection

- ✓ Fast Results in 10 minutes
- ✓ Simple Procedure
- ✓ Suitable for all Geographies
- ✓ High Accuracy

Analytes	Sensitivity	Specificity	Accuracy
P.f.	>99.9%	>99.9%	>99.9%
P.v.	98.2	>99.9%	99%
Pan	>99.9%	>99.9%	>99.9%



- Malaria is caused by Plasmodium sp. (parasite) which is spread by mosquitos.
- Malaria is widespread in the tropical and subtropical regions that exist in a broad band around the equator. This includes much of Sub-Saharan Africa, Asia, and Latin America.
- Symptoms usually begin ten to fifteen days after being bitten by an infected mosquito. If not properly treated, people may have recurrences of the disease months later.
- CITEST Malaria Rapid Tests are rapid chromatographic immunoassay for the qualitative detection of P.f.-specific HRP - II, P.v.-specific LDH and Pan-malarial Aldolase (P.f., P.v., P.o. and P.m.) antigens in whole blood with accurate results.

#### Convenient Procedure



#### Ordering Information

Cat. No.	Product	Specimen	Format	Pack	CE Status
IMPVF -402	Malaria P.f./P.v. /Pan Rapid Test Cassette	WB	Cassette	25T	CE