




CPV+CDV Antigen Combo VETERINARY TEST



VIDP-525

PRODUCT INFORMATION

Canine Parvovirus (CPV)+Canine Distemper Virus (CDV) Antigen Combo

Specimen: Secretions/Feces/Vomit
Format: Cassette Result Time: 5-10 min
Pack: 10T Certificate: /

	Sensitivity	Specificity	Accuracy
CDV Ag	97.96%	97.50%	97.75%
CPV Ag	95.74%	>99.99%	97.56%

Canine distemper is a viral disease that affects a wide variety of mammal families. The disease is highly contagious via inhalation. In domestic dogs, while the acute generalized form of distemper has a high mortality rate.

Canine parvovirus is a contagious virus mainly affecting dogs. CPV is highly contagious and is spread from dog to dog by direct or indirect contact with their feces. Once infected, the white blood cell level falls, further weakening the dog. Younger animals have worse survival rates.

Phone: 012-313 3081

Email: sales@vigagemilang.com.my

B-G-8, RAWANG IDAMAN

JALAN IDAMAN 1, TAMAN RAWANG IDAMAN,
48000 RAWANG, SELANGOR, MALAYSIA