



# FIV Ab + FELV Ag Combo VETERINARY TEST



**VIIE-CT425**

## PRODUCT INFORMATION

Feline Immunodeficiency Virus (FIV) Antibody & Feline Leukemia Virus (FeLV) Antigen Combo

Specimen: WB/S/P

Format: Cassette

Pack: 10 T

Result Time: 5-10min

FIV Ab:

Sensitivity: 97.50%

Specificity: 96.67%

Accuracy: 96.88%

FELV Ag:

Sensitivity: 92.86%

Specificity: >99.99%

Accuracy: 97.06%

Feline immunodeficiency virus (FIV) is a Lentivirus that affects cats worldwide. FIV progresses through similar stages to HIV in humans. The initial stage is fairly short and is followed by the asymptomatic stage. Finally, the cat progresses into the final stage (AIDS stage), where in the cat is extremely susceptible to secondary diseases that inevitably are the cause of death.

Feline leukemia virus (FeLV) is a retrovirus that infects cats. If not defeated by the animal's immune system, the virus weakens the cat's immune system which can lead to diseases which can be lethal.

Phone: 012-313 3081

Email: [sales@vigagemilang.com.my](mailto:sales@vigagemilang.com.my)

B-G-8, RAWANG IDAMAN

JALAN IDAMAN 1, TAMAN RAWANG IDAMAN,

48000 RAWANG, SELANGOR, MALAYSIA