


SWAB SUMMARY

COLLECTION
TRANSPORT
RECOVERY



BETTER DIAGNOSTICS BEGIN WITH BETTER SAMPLES

And better samples are obtained by using proper sample collection methods. Through market research and clinical evaluations, Babio has designed and developed flocked swab that optimize the efficiency of the target analyte collection with improved patient comfort. Babio Swabs head is made of artificial fiber material, and act like a soft brush and allow improved collection of cell samples. It cannot disperse and entrap the precious sample like traditional fiber wound swabs. In stark contrast, the entire sample stays close to the surface for fast and complete elution. The molded breakpoint allows safe and easy breaking of the swab stick. Several breakpoint distance options available for fitting different tubes. Babio Swab use short fiber strand which create a thin absorbent layer. The entire sample stays close to the surface for fast and complete elution.



SWAB ADVANTAGE







- 1 ERGONOMIC AND ANATOMIC DESIGN
—improves patient comfort and specimen collection.
- 2 RAPID AUTOMATIC ELUTION
—for instant and spontaneous release of sample into liquid media.
- 3 MEASURABLE AND CONSISTENT
— uptake and transfer from patient for a quantitative measurement and improved test sensitivity
- 4 PATIENT FRIENDLY SWAB
— have anatomical designs and a soft brush like texture for rapidly and efficiently dislodging cells.
- 5 SWAB ARE COMPATIBLE WITH MULTIPLE APPLICATIONS AND PLATFORMS
—such as rapid antigen testing, EIA, molecular-based assays, DFA, cytology testing, bacteriology and virology culture.
- 6 SAFE AND CONVENIENT MOLDED BREAKPOINT.



SWAB STERILE

Product Description						
Break Point Distance (From Tip End)	80mm	80mm	80mm	80mm	80mm	80mm
Product Code	N15080	NT15080	NTP15080	T15080	TT15080	TTP15080
Product Description						
Break Point Distance (From Tip End)	30mm	30mm	30mm	No Break Point	No Break Point	No Break Point
Product Code	T15030	TT15030	TTP15030	M15000	MT15000	MTP15000

SWAB APPLICATION

		
Use Alone	Use With VTM	Use With Virus Transport Ki (tInactivated)
		
Use With Media	Use With Liquid Amies Media	Use With Liquid Other Media